




HA(Δ TM)(H₁N₁/A/New Caledonia/20/99)

Product Name	HA(Δ TM)(H ₁ N ₁ /A/New Caledonia/20/99)
Catalog No	IT-003-001 Δ TMp
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(Δ TM) (H ₁ N ₁ /A/New Caledonia/20/99) (a.a.18-530)(Genebank No. AAP34324)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM)(H₁N₁/A/New Caledonia/20/99)

Note: