


HA1(H1N2)(A/swine/Gent/7625/1999)



Product Name	HA1(H1N2)(A/swine/Gent/7625/1999)
Catalog No	IT-003-001SW2p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA1(H1N2)(A/swine/Gent/7625/1999) (a.a.18-336) (Genebank Accession No. AAT00437)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA1(H1N2) (A/swine/Gent/7625/1999)