



HA(Δ TM)(H2N2)(A/Canada/720/05)

Product Name	HA(Δ TM)(H2N2)(A/Canada/720/05)
Catalog No	IT-003-0031p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA (H2N2)(A/Canada/720/05) protein (a.a.13-526) (Genebank No.AAY28987).
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM)(H2N2)(A/Canada/720/05) protein

Note: