

H3(Δ TM)(H3N2)(A/Guizhou/54/89)



Product Name	H3(Δ TM)(H3N2)(A/Guizhou/54/89)
Catalog No	IT-003-00412 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged H3(Δ TM)(H3N2)(A/Guizhou/54/89) protein (a.a.17-529) (Genebank #:BAA08716)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: H3(Δ TM)(H3N2)(A/Guizhou/54/89)