




HA1(H3N2)(A/Christchurch/4/85)

Product Name	HA1(H3N2)(A/Christchurch/4/85)
Catalog No	IT-003-00414p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA1(H3N2)(A/Christchurch/4/85)(1-329aa; GenBank# AAB69838)
Applications	WB standard, antibody ELISA, etc
Purity	≥ 95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA1(H3N2)(A/Christchurch/4/85)

Note: