


H3(Δ TM)(H3N2) (A/Mississippi/1/85)



Product Name	H3(Δ TM)(H3N2)(A/Mississippi/1/85)
Catalog No	IT-003-00415 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged H3(Δ TM)(H3N2)(A/Mississippi/1/85) protein (a.a.17-529) (Genebank #: AAA62338)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: H3(Δ TM)(H3N2)(A/Mississippi/1/85)