


H3(Δ TM)(H3N2)(A/Victoria/3/75)



| | |
|----------------|--|
| Product Name | H3(Δ TM)(H3N2)(A/Victoria/3/75) |
| Catalog No | IT-003-00417 Δ TMp |
| Source | Recombinant protein purified from 293 cell culture |
| Formulation | Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS |
| Protein | C-terminal 6xHis tagged H3(Δ TM)(H3N2)(A/Victoria/3/75) protein (a.a.17-530) (Genebank #: P03435) |
| Applications | WB standard, antibody ELISA, immunogen, etc |
| Purity | >95% |
| Storage | Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |



E: H3(Δ TM)(H3N2)(A/Victoria/3/75)