

H3(Δ TM)(H3N2) (A/Hong Kong/8/68)

Product Name	H3(Δ TM)(H3N2)(A/Hong Kong/8/68)
Catalog No	IT-003-00418 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged H3(Δ TM)(H3N2)(A/Hong Kong/8/68) protein (a.a.17-529) (Genebank #: AAK51718)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: H3(Δ TM)(H3N2)(A/Hong Kong/8/68)