


# H3( $\Delta$ TM)(H3N2) (A/England/42/72)



---

Product Name	H3( $\Delta$ TM)(H3N2)(A/England/42/72)
Catalog No	IT-003-00420 $\Delta$ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 $\mu$ g of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged H3( $\Delta$ TM)(H3N2)(A/England/42/72) protein (a.a.17-530) (Genebank #: ABQ58930)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: H3( $\Delta$ TM)(H3N2)(A/England/42/72)