

HA₁(H₃N₂)(A/VICTORIA/210/2009)



Product Name	HA ₁ (H ₃ N ₂)(A/VICTORIA/210/2009)
Catalog No	IT-003-00421p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA ₁ (H ₃ N ₂)(A/VICTORIA/210/2009) protein(a.a. 19-350)(Genebank #: EPI272062)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA₁(H₃N₂)(A/VICTORIA/210/2009)