


HA(Δ TM)(H3N2) (A/Switzerland/9715293/2013)



Product Name	HA(Δ TM)(H3N2)(A/Switzerland/9715293/2013)
Catalog No	IT-003-00428 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA(Δ TM)(H3N2)(A/Switzerland/9715293/2013) (see below for protein sequence)
Applications	WB standard, antibody ELISA, etc
Purity	$\geq 95\%$
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Protein SEQ:

qklpgndstatlclghhavpngtivrktitndrievtnatelvqnssigeicdsphqildgenctlidallgdpqcdgfnkkw
dlfverskaysncypdyvpdyaslrlsvassgtlefnnesfnwagvtqngtssscirgnsnsffsrlnwlthlnskypalnvmt
pnneqfdklyiwgvhhpgtdkdqiflyaqssgritvstkrsqqavipnigrprirdipsrisiywtivkpgdillinstgnliap
rgyfkirsgkssimrdsapigkcksecitpngsipndkpfqvnvritygacpryvkqstlklatgmrnvperqtrgifgaiagf
iengwegmvdgwygfrhqnsgrgqaadlkstqaaidqingklnrligktnekfhqiekfsevegriqdlekyvedtkid
lwsynaellvalenqhtidldtsemnklfektkkqlrenaedmngcfcfkiyhkcdnacigsirnngtydhdvyrdealnnrfq
ikgvelksgykd