



# HA<sub>1</sub>(H<sub>3</sub>N<sub>2</sub>)(Influenza A)

---

Product Name	HA <sub>1</sub> (H <sub>3</sub> N <sub>2</sub> ) (A/Brisbane/10/2007)
Catalog No	IT-003-0042p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA <sub>1</sub> (H <sub>3</sub> N <sub>2</sub> ) (A/Brisbane/10/2007) protein (a.a.17-345) (Genebank No. EU199248)
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing.</u> The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

