

# HA1(A/Aichi/2/1968(H<sub>3</sub>N<sub>2</sub>))



---

Product Name	HA1(A/Aichi/2/1968(H <sub>3</sub> N <sub>2</sub> ))
Catalog No	IT-003-0043p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA1(A/Aichi/2/1968(H <sub>3</sub> N <sub>2</sub> )) protein (a.a. 17-346)(Genebank No. AAA43178).
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing.</u> The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

