




HA(Δ TM)(H3N2/Wisconsin/15/2009)

Product Name	HA(Δ TM) (H3N2/Wisconsin/15/2009)
Catalog No	IT-003-0044 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(Δ TM) (H3N2/Wisconsin/15/2009) protein (a.a.18-529) (Genebank #: ACR32998)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM) (H3N2/Wisconsin/15/2009)