



## HA1(H3N2)(A/Panama/07/99)

---

Product Name	HA1(H3N2)(A/Panama/07/99)
Catalog No	IT-003-0047p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA1(H3N2)(A/Panama/07/99)(1-329aa; GenBank# ABH01004)
Applications	WB standard, antibody ELISA, etc
Purity	≥ 95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

