

# HA<sub>1</sub>(A/Pennsylvania/14/2010(H3N2))



---

Product Name	HA <sub>1</sub> (A/Pennsylvania/14/2010(H3N2))
Catalog No	IT-003-004SW4p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA <sub>1</sub> (A/Pennsylvania/14/2010 (H3N2)) (a.a.17-345)(Gene Bank Accession No.ADD21560) (human isolate of SOIV H3N2)
Applications	Immunogen, antibody ELISA, standard, etc.
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA<sub>1</sub>(A/Pennsylvania/14/2010(H3N2))