

HA₁(A/Ontario/RV1273/2005(H3N2))



Product Name	HA ₁ (A/Ontario/RV1273/2005(H3N2))
Catalog No	IT-003-004SW5p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA ₁ (A/Ontario/RV1273/2005(H3N2))(a.a.17-345)(Gene Bank Accession No. ABF17954) (human isolate of SOIV H3N2)
Applications	Immunogen, antibody ELISA, standard, etc.
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA₁(A/Ontario/RV1273/2005(H3N2))