

# HA1(H3N2v)(A/Indiana/07/2012)



---

Product Name	HA1(H3N2v) (A/Indiana/07/2012)
Catalog No	IT-003-004SW6p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA1(H3N2v) (A/Indiana/07/2012)(a.a.13-342)(Genbank Accession No. ACN98096) ,which is a strain recently transmitted to human.
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

