

HA1(H3N2)(A/swine/Flanders/1/1998)



Product Name	HA1(H3N2)(A/swine/Flanders/1/1998)
Catalog No	IT-003-004SW7p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA1(H3N2)(A/swine/Flanders/1/1998) (a.a.13-342)(Genebank Accession No. ACN98096)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

