




H₅ (H₅N₁)(A/Vietnam/1203/2004)

Product Name	H ₅ (H ₅ N ₁)(A/Vietnam/1203/2004)
Catalog No	IT-003-0051p-PBS
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged H ₅ (H ₅ N ₁)(A/Vietnam/1203/2004) protein (a.a.18-530; Genebank# AAW80717) with a deletion at its cleavage site (ΔRRRKKR). Majority of the protein exists as trimer/oligomer.
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: purified H₅ (H₅N₁)(A/Vietnam/1203/2004)