


HA₁(H₅N₁)(A/Egypt/N03072/2010)



Product Name	HA ₁ (H ₅ N ₁)(A/Egypt/N03072/2010)
Catalog No	IT-003-00520p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA ₁ (H ₅ N ₁)(A/Egypt/N03072/2010) protein(a.a. 18-345)(Genebank #: EPI255379)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA₁(H₅N₁)(A/Egypt/N03072/2010)