


# HA1 (A/chicken/Bangladesh/11rs1984-30/2011) (H5N1)



Product Name	HA1(A/chicken/Bangladesh/11rs1984-30/2011)(H5N1)
Catalog No	IT-003-00525p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA1(A/chicken/Bangladesh/11rs1984-30/2011) (H5N1) protein (a.a. 20-338) (see amino acid sequence below)(Genebank Accession No. AEQ50043)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



## Protein SEQ

nnsteqvdtimkntvtvthaqdilekthngklcdldgvpkplilrdcsvagwllgnpmcdefinvpewsyivekanpandl  
cypgnlneyeelkhllsrinhfekiipksswsdheaslgvssacpylgtpsffrnvwlikknntyptikrsyntnqedlli  
lwgihhpneaeqiqllyqnpttyisvgtstlnqrlvpkiatrskvngqsgmrdffwtilkpddainfesngnfiapeyaykiv  
kkgdsaimksdveygncntkccqtpvgainssmpfhnihltigecpkyvksnklvlatglrnsalrekrk