


HA1(A/turkey/Germany-MV/R2472/2014)(H5N8)



Product Name	HA1(A/turkey/Germany-MV/R2472/2014)(H5N8)
Catalog No	IT-003-00528p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA1(A/turkey/Germany-MV/R2472/2014)(H5N8) (see below for protein sequence)
Applications	WB standard, antibody ELISA, etc
Purity	≥ 95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA1(A/turkey/Germany-MV/R2472/2014)(H5N8)

Note:

Protein SEQ:

dqicigyhannstkqvdtimeknvtvthaqdilekthngklcdlngvklilkdcsvagwllgnpmcdefirvpewsyive
ranpandlcypgtlndyeelkhllsrinhfektliipksswvnhetslgvsaacpyqgassffrnvwwlikkndayptikisyn
ntredllilwgihhpnaaeqtnlyknpdtyvsvgtstlnqrlvpkiatrsqvngqrgrmdffwtilkpndaihfesngnfia
peyaykivkkgdstimkseveyghcntkcqtpigainssmpfhnihpligtgecpkyvksnklvlatglmsplreg