




# HA( $\Delta$ TM)(H<sub>5</sub>N<sub>1</sub>)

---

Product Name	HA( $\Delta$ TM)(A/Indonesia/5/2005(H5N1))
Catalog No	IT-003-0052 $\Delta$ TMp
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 $\mu$ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA( $\Delta$ TM)(H <sub>5</sub> N <sub>1</sub> ) (A/Indonesia/5/2005) protein (a.a.17-530 ) (Genebank No. <a href="#">ABP51969</a> )
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA( $\Delta$ TM)(H<sub>5</sub>N<sub>1</sub>)(A/Indonesia/5/2005) protein

Note: