




# H<sub>5</sub> (H<sub>5</sub>N<sub>1</sub>)(A/Indonesia/5/2005)

---

Product Name	H <sub>5</sub> (H <sub>5</sub> N <sub>1</sub> )(A/Indonesia/5/2005)
Catalog No	IT-003-0052p-PBS
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged H <sub>5</sub> (H <sub>5</sub> N <sub>1</sub> )(A/Indonesia/5/2005) protein (a.a.18-530; Genebank# ABP51969) with a deletion at its cleavage site (ΔRRRKKR). Majority of the protein exists as trimer/oligomer.
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: purified H<sub>5</sub> (H<sub>5</sub>N<sub>1</sub>)(A/Indonesia/5/2005).