


HA1(A/northern pintail/Washington/40964/2014)(H5N2)



Product Name	HA1(A/northern pintail/Washington/40964/2014)(H5N2)
Catalog No	IT-003-00531p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA1(A/northern pintail/Washington/40964/2014)(H5N2)(see below for its protein seq)
Applications	WB standard, antibody ELISA, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Protein SEQ:

dqicvgyhannstervdtimeknvtvthaqdilekthngklcnlmgvklpledcsvagwllgnpmcdkflnvsewsyive
kagpanglcypgdfndyeelkhllsrinhfekieipknhwsnhntsgvssacsylenpsffrnvvltkknntyptikvny
tnanqkdllvlwgihhpnaeqkmiyqnsntyvsvgtstlnqrlvpkiatrskvngqsgridffwtilkpndtinfdsngnf
iapkyaykivkkgdsaimkselkydncntkcqtpigainssmpfhnihpltigecpkyvksnrlvlatglmapqiegrr