



HA₁(H₅N₁)(A/Hong Kong/156/97)

Product Name	HA ₁ (H ₅ N ₁) (A/Hong Kong/156/97)
Catalog No	IT-003-0056p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA ₁ (H ₅ N ₁) (A/Hong Kong/156/97) protein (a.a.18-346) (Genebank No. AAC32088)
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA₁(H₅N₁) protein

Note: