

HA2(A/Vietnam/1203/2004)(H5N1)



Product Name	HA2(A/Vietnam/1203/2004)(H5N1)
Catalog No	IT-003-0058p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA2 (A/Vietnam/1203/2004)(H5N1) (a.a.366-531) (Genebank No. AY651334).
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

