



# HA<sub>1</sub>(H<sub>4</sub>N<sub>6</sub>)(Influenza A)

---

|                |  |
|----------------|--|
| Product Name   | HA <sub>1</sub> (H <sub>4</sub> N <sub>6</sub> )   |
| Catalog No     | IT-003-006p  |
| Source         | Viral protein purified from 293 cell culture   |
| Formulation    | Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.  |
| Protein        | 6xHis tagged HA <sub>1</sub> (H <sub>4</sub> N <sub>6</sub> ) (A/mallard/Ohio/657/2002) protein (a.a.17-343) (Genebank No. ABI47995) |
| Applications   | WB standard, antibody ELISA, etc   |
| Purity         | >95%   |
| Storage        | Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.                   |
| Use Limitation | For research use only, not for use in diagnostic procedures.   |



E: HA<sub>1</sub>(H<sub>4</sub>N<sub>6</sub>) protein

Note: