


# HA1 (A/chicken/Jalisco/CPA1/2012) (H7N3)



Product Name	HA1 (A/chicken/Jalisco/CPA1/2012(H7N3))
Catalog No	IT-003-0071p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA1 (A/chicken/Jalisco/CPA1/2012) (H7N3)(a.a.19-344) (See its amino acid sequence below; Genebank #: AFN85519)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



## Protein SEQ:

dkiclgghavangtkvntltergievvnatetvetantkkictqgkrptdlgqcgllgtligppqcdqfledadliierregtdvc  
ypgkftneeslrqilrgsggidkesmgftysgirtngatsacrrsgssfyamkwllnsdnaafpqmtksyrnprnkpalii  
wgvhhsqsateqtklysgnklitvgsskyqqsfpspgarpqvnngqsgridfhwllldpndvtftfngafiapdrasffrge  
sigvqsdvpldsgegdcfhnggtivsslpfqinprtvgkcpnyvkqtslllatgmrnvpnpkdrksrhr