


# HA1(H7N9)(A/Shanghai/1/2013)



---

Product Name	HA1(H7N9)(A/Shanghai/1/2013)
Catalog No	IT-003-0073p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA1 (H7N9)(A/Shanghai/1/2013)(a.a.19-339) (GISAID #: EPI439486). This is the HA1 protein from the outbreak of the novel H7N9 virus in China since late March, 2013.
Applications	Protein standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA1(H7N9)(A/Shanghai/1/2013)