

HA1 (A/Hong Kong/1073/99) (H9N2)



Product Name	HA1(A/Hong Kong/1073/99)(H9N2)
Catalog No	IT-003-0091p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA1(A/Hong Kong/1073/99)(H9N2) protein (a.a.19-339) (Genbank #: NP_859037)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

