

# HA<sub>1</sub>(H<sub>9</sub>N<sub>2</sub>)( A/Hong Kong/33982/2009)



---

Product Name	HA <sub>1</sub> (H <sub>9</sub> N <sub>2</sub> )( A/Hong Kong/33982/2009)
Catalog No	IT-003-0094p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA <sub>1</sub> (H <sub>9</sub> N <sub>2</sub> )( A/Hong Kong/33982/2009) protein(a.a. 17-339)(Genebank #: EPI241740)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA<sub>1</sub>(H<sub>9</sub>N<sub>2</sub>)( A/Hong Kong/33982/2009)