


HA1 (A/chicken/India/IVRI-0011/2011) (H9N2)



Product Name	HA1(A/chicken/India/IVRI-0011/2011)(H9N2)
Catalog No	IT-003-0095p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA1(A/chicken/India/IVRI-0011/2011)(H9N2)protein (a.a. 16-335) (see amino acid sequence below)(Genebank Accession No. AEP84155)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Protein SEQ

dkicighqstnstetvdtltetnvpvthakellhtehngmlcatnlqgpilntctieglygnpscdlllgrewsyiverpsav
ngtcypgnvenleelrtlfssssyriqifpdtiwnvftgtkskscdsfyrnmrwlqknglypvqdaqytnnqgrdlfvw
gihhpptdaqtnlytrtdtttsvttenldrtfkpligrprplvngligrinyysvltkpgqtlrvrsngnliapwyghvlsghg
rilktdlksngnevvcqctekgglnstlpfhniskyafgtcpsyigvksklaiqlrnvparssr