

HA1(A/spot-billed duck/Jiangxi/33/2011) (H9N2)



Product Name	HA1(A/spot-billed duck/Jiangxi/33/2011)(H9N2)
Catalog No	IT-003-0096p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA(A/spot-billed duck/Jiangxi/33/2011)(H9N2) (a.a.19-339) (GISAID# EPI439646)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA1(A/spot-billed duck/Jiangxi/33/2011) (H9N2)

Protein SEQ:

dkicighqstnstetvdtltetnvpvthakellqtehnmgmlcatnlghplildtctiegliygnpdscdllegrewsyiverpsav
ngtcypgnvenleelrtlfssassyqriqifpdttnvntytgtsrscsgsfyrsmrwlqksgfyypvdakyttnrgknilfvw
gihhphtytdqttlirndtttsvateelnrifkpvivprllvngqqgridyywsvlkgqtlrvsnngniapwyghvlsggsh
grilktksksgncvvcqtekggfntlpfhniskyafgtcpsykvnsklaiqlrnparsnrg