

HA1(A/duck/Guangdong/E1/2012)(H10N8)



Product Name	HA1(A/duck/Guangdong/E1/2012)(H10N8)
Catalog No	IT-003-0111p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA1(A/duck/Guangdong/E1/2012)(H10N8) protein (a.a.18-341; Genbank# JQ924786).
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

