

HA1(A/quail/Hong Kong/1721-30/99(H6N1))



Product Name	HA1(H6N1)(A/quail/Hong Kong/1721-30/99)
Catalog No	IT-003-0211p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA1(H6N1)(A/quail/Hong Kong/1721-30/99) protein (a.a.17-346) (Genbank #: CAC84981)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA1(H6N1)(A/quail/Hong Kong/1721-30/99)