


H17(Δ TM)(A/little yellow-shoulderedbat /Guatemala/153/2009(H17N10))



Product Name	H17(Δ TM)(A/little yellow-shoulderedbat/Guatemala/153/2009(H17N10))
Catalog No	IT-003-022 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged H17(Δ TM)(A/little yellow-shoulderedbat /Guatemala/153/2009(H17N10)) (Genbank Accession No. AFC35418.1)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: H17(Δ TM)(A/little yellow shoulderedbat /Guatemala/153/2009(H17N10))