


M1(A/Perth/16/2009)(H3N2)



Product Name	M1(A/Perth/16/2009)(H3N2)
Catalog No	IT-003-039Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein in PBS with 8M Urea. Reconstitute the protein with 100 µl of Millipore water.
Protein	C-terminal 6xHis tagged M1(A/Perth/16/2009)(H3N2) protein (a.a. 1-252) (See below for sequence).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Before reconstitution, stable for 1 year at -20°C from the date of shipment. After reconstitution, stable for a month at 4°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Protein SEQ:

MSLLTEVETYVLSIVPSGPLKAEIAQRLEDVFAGKNTDLEALMEWLKTRPILSPLT
KGILGFVFTLTVPSERGLQRRRFVQNALNGNGDPNNMDKAVKLYRCLKREITFH
GAKEIALSYSAGALASCMGLIYNRMGAVTTEVAFGLVCATCEQIADSQHRSHRQ
MVATTNPLIKHENRMVLA STTAKAMEQ MAGSSEQAAEAMEIASQARQM VQAM
RAIGTHPSSSTGLRDDLLENLQTYQKRMGVQMQRFK