


M1 (A/California/07/2009)(H1N1)



Product Name	M1 (A/California/07/2009)(H1N1)
Catalog No	IT-003-040Ep
Source	Viral protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein (1 mg/ml) in PBS. Reconstitute the protein with 100µl of Millipore water
Protein	C-terminal 6xHis tagged M1 (A/California/07/2009)(H1N1) matrix protein (Genbank accession no. AFM72833)
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Store at -80°C; <u>AVOID FREEZING/THAWING</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: M1 (A/California/07/2009)(H1N1) protein

Note: