

# M1 (B/Malaysia/2506/2004)



---

Product Name	M1 (B/Malaysia/2506/2004)
Catalog No	IT-003-042Ep
Source	Recombinant protein expressed in E.coli
Formulation	Each vial contains 100 µg of lyophilized protein in PBS with 8M Urea. Reconstitute the protein with 100 µl of Millipore water
Protein	C-terminal 6xHis tagged M1 (B/Malaysia/2506/2004) (a.a.1-248)(Genebank Accession No. ACR15733)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at -80°C. Stable for 1 month after reconstitution at 4°C. For long term storage after reconstitution, aliquot and freeze at -80°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

