


HA1(B/YAMANASHI/4/2012)



Product Name	HA1(B/YAMANASHI/4/2012)
Catalog No	IT-003-B10p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA1(B/YAMANASHI/4/2012) protein (a.a. 16-362) (see amino acid sequence below)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Protein SEQ

driectgitsnsphvvtatqgevnvtgviplttptksyfanlkgtrtrgkclpcdlnctldvalgrpmcvgtppsakasilhe
vrpvtsgcfpimhdrtkirqlpnllrgyekirlstqnvidaekapggpyrlgtsgscpnatskigffatmawavpkdnyknt
npltevpyictggedqitvwgfhdsnktqmkslygdsnpqkftssangvtthyvsqjgdfpdqtedgglpqsgrivvdy
mmqkpgktgativyqrgvllpqkvwcasgrskvikgslpligeadclhekygglnskpytgehakaigncpiwvktpl
klangkyrppakllkerg