

# HA( $\Delta$ TM)(B/Harbin/7/1994)



Product Name	HA( $\Delta$ TM)(B/Harbin/7/1994)
Catalog No	IT-003-B15 $\Delta$ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 $\mu$ g of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA( $\Delta$ TM)(B/Harbin/7/1994) protein (a.a. 30-547) (Genebank Accession No. ACR15721)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

