

# HA1(B/Bangladesh/5945/2009)



---

|                |  |
|----------------|--|
| Product Name   | HA1(B/Bangladesh/5945/2009)  |
| Catalog No     | IT-003-B17p  |
| Source         | Recombinant protein purified from 293 cell culture   |
| Formulation    | Each vial contains 100 µg of purified protein (1 mg/ml) in PBS   |
| Protein        | C-terminal 6xHis tagged HA1(B/Bangladesh/5945/2009) (a.a.18-362) (Genbank #: AFH58447).                              |
| Applications   | Protein standard, antibody ELISA, immunogen, etc   |
| Purity         | >95%   |
| Storage        | Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures.   |

