



HA₁(A/Swine/Wisconsin/136/97 (H₁N₁))

Product Name	HA ₁ (A/Swine/Wisconsin/136/97(H ₁ N ₁))
Catalog No	IT-003-SW6p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA ₁ (A/Swine/Wisconsin/136/97(H ₁ N ₁)) protein (a.a.18-342) (Genebank No. AAF87275)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA₁(A/Swine/Wisconsin/136/97(H₁N₁))