



# Nef (HIV-1/Clade B)

---

Product Name	Nef (HIV-1/Clade B)
Catalog No	IT-001-014Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of protein in PBS with 8M Urea.
Protein	N-terminal GST- and C-terminal 6xHis tagged HIV-1 Nef (clade B) protein (GenBank No. <a href="#">K03455</a> )
Applications	WB standard, antibody ELISA, etc
Purity	>90%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

