



# Nuclear antigen-1 (EBV)

---

|                |  |
|----------------|--|
| Product Name   | Nuclear antigen-1 (EBV)  |
| Catalog No     | IT-005-051Ep   |
| Source         | Recombinant protein purified from E.coli   |
| Formulation    | Each vial contains 100 µg of protein in PBS with 8M Urea.  |
| Protein        | N-terminal GST and C-terminal 6xHis tagged mosaic protein containing the HHV-4 EBNA regions 1-90 and 408-498 amino acids. (GenBank Accession No: AGZ95187).  |
| Applications   | WB standard, antibody ELISA, immunogen, etc  |
| Purity         | >85% , based on SDS PAGE   |
| Storage        | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures.   |

