


# NS5a (HCV/Genotype-1a/Con)



Product Name	NS5a (HCV/Genotype-1a/Con)
Catalog No	IT-004-078Ep
Source	Viral protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein (1 mg/ml) in PBS.
Protein	GST-fused HCV NS5a Genotype-1a consensus immunodominant regions (aa 2213-2332).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>90% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Protein SEQ:

katctanhdsdpaieianllwrqemggnitrvesenkvvildsfdplvaeedereisvpaeilrksrrfapalpiwarpdynp  
pilletwkkpdyepvvhgcplpppqspvpvpprkr