


# p24 (HIV-1/Clade B)(HXBc2)



Product Name	p24 (HIV-1/Clade B)(HXBc2)
Catalog No	IT-001-017Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein in PBS with 8M Urea.
Protein	6xHis tagged p24 (HIV-1/Clade B)(HXBc2) antigen (a.a.134-363) (GenBank Accession No. AAB50258).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: p24 (HIV-1/Clade B)(HXBc2)