

HA(Δ TM)(A/Bolivia/559/2013)(H1N1)



Product Name	HA(Δ TM)(A/Bolivia/559/2013)(H1N1)
Catalog No	IT-003-00108 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA(A/Bolivia/559/2013)(H1N1)(a.a.18-529)(GISAID No. EPI466837).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	\geq 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

